

WELCOME TO RESPIMAT®

Respimat® has been designed for effortless inhalation. At each prescription, you can help your patients prepare Respimat® so they only need to focus on daily use.¹

Also, Respimat® Instructional Tear Sheets have been provided for you to give to patients.

For Respimat® placebo devices, please contact your local Boehringer Ingelheim representative.*

Thank you.

*Placebo Respimat® is shown for demonstration purposes only.

*With SPIRIVA® Respimat®, the medication is administered as two puffs once daily.



RESPIMAT®
Unique mist, reliable delivery.

RESPIMAT® IN DAILY USE: JUST REMEMBER TOP²



TURN

- Keep the cap closed.
- **TURN** the clear base in direction of the arrows on label until it clicks.



OPEN

- **OPEN** the cap until it snaps fully open.



PRESS

- Breathe out slowly and fully.
- Close your lips around the mouthpiece.
- While taking a slow, deep breath through your mouth, **PRESS** the dose release button and continue to breathe in.
- Hold your breath as long as comfortable.
- Close the cap.

¹With SPIRIVA® Respimat®, the medication is administered as two puffs once daily.



RESPIMAT®
Unique mist, reliable delivery.

RESPIMAT® PREPARATION: ONLY BEFORE FIRST USE²



- Keep the cap closed.
- Press the safety catch and pull off the clear base.



- Keep the cap closed.
- Insert narrow end of cartridge into inhaler.
- Place inhaler on a firm surface and push down firmly until it clicks into place.



- Keep the cap closed.
- Replace the clear base.



- Keep the cap closed.
- Turn clear base in direction of the arrows on label until it clicks.



- Open the cap until it snaps fully open.



- Point Respimat® toward ground and press the dose release button.
 - Close the cap.
- Repeat steps 4–6 until a cloud is visible. After a cloud is visible repeat steps 4–6 three more times.

For a video tutorial, please visit <http://interactive.respimat.com> or www.respimat.com.
For more details please refer to the package leaflet of your Respimat® product.

References: 1. Ciciliani A et al. *Respiratory Drug Delivery*. 2015;2:357-362. 2. Placebo Respimat® inhaler package leaflet; instructions for demonstration. Boehringer Ingelheim International GmbH. July 2016.



© Boehringer Ingelheim International GmbH. All rights reserved. August 2016

RESPIMAT®
Unique mist, reliable delivery.